



Devices to help apply stockings

The use of application aids not only assist with getting elastic compression stockings on, they can help extend the effectiveness of the garments.

Take advantage of these devices to help apply your stockings today.



medi textile gloves
Textured to allow for easy fabric control.



big butler
For the patient with a full calf.

How to care for your stockings

Proper care ensures not only that your stockings last, but also that they continue to provide maximum therapeutic benefits in your vein disease treatment.

The below medi products can assist you in proper care of your garments. Remember that clean compression equals effective compression.



medi night care cream for nighttime use specifically formulated to provide the ingredients your skin needs to recuperate and regenerate overnight. Simply apply this soothing creme to your legs after removing your compression stockings.



it stays This odorless, hypoallergenic and **easy-to-wash adhesive** securely and comfortably holds knee-highs, thigh-highs and arm sleeves in place.



medi care™ wash Extend the life of your compression therapy products by using **medi care™ wash** every day. Our pure, organic, and biodegradable solution thoroughly removes the daily accumulation of perspiration and body oils from compression garments.

medi USA
6481 Franz Warner Parkway
Whitsett, NC 27377
T 1.800.633.6334
F 888-570-4554
info@mediusa.com
mediusa.com

PO150 REV G

medi



Compression Garment Care
Everything you need to know to use compression

medi. I feel better.

WE ARE HERE TO HELP!

Let medi® be your personal Compression Therapy Coach.

Join the team and receive compression tips, special offers, expert advice and more to help you meet your therapy goals!




Sign up NOW!



Your doctor has prescribed compression stockings

to help alleviate symptoms of your specific condition. To ensure the best therapeutic results, medical stockings should be worn every day and replaced every six months. Wear exactly as prescribed. Follow your doctor's directions for when and how long to wear your stockings.

Use your smart phone to scan the QR codes to the right!




How to apply stockings

These techniques will help you apply your stockings easily and correctly. Keep in mind that swelling is less pronounced in the morning, so applying your mediven®, duomed® & Rejuva® stockings as soon as you get out of bed will make application easier.

Standard techniques

Applying closed-toe stockings.



Video available for closed-toe stockings application



Applying Open-toe Stockings

Applying open-toe stockings using an easy-on slipper.



Video available for open-toe application



Check out these other video tutorials for application and removal devices



butler
For use with all knee-high and thigh-high products.



Application with butler device



medi butler off
Enables easy removal by extending your reach.



Removal with butler device



medi 2in1
Aid for open and closed toe garments.



Removal of 2-in-1